

eka-Lpt cushions with cover of incontinence material

CLEANING

With wipe disinfection:

To further ensure the product properties of the material we recommend the following products for wipe disinfection:

- Dimozon pur 4%
- Formaldehyd 30%
- Incidin Plus 2%
- Acryl-des 100
- Descogen 10%
- Mikrozid AF liquide 100%
- Perform 6%
- TPH protect 4%

Proceed as follows for wipe disinfection:

Spray the product with the wipe disinfectant or apply with a cloth and allow to dry.

Please observe the manufacturer's instructions!

For further cleaning instructions, please refer to the label on the product.